Arthrex ACP® Double Syringe

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Ordering Information

<i>Item number</i>	<i>Description</i>
ABS-10011	Arthrex ACP Kit, Series I
ABS-10012	Arthrex ACP Kit, Series II
ABS-10014	Arthrex ACP Double Syringe
ABS-10020	Centrifuge Hettich Rotofix 32 with Swing Out Rotor, 220 V
ABS-10022	Bucket with Screw Cap for Centrifuge (spare)
ABS-10023	Screw Cap for ABS-10022 (spare)
ABS-10027	Counterweight for Centrifugation of ACP-Double Syringe, 15 ml
ABS-10050	ViscoGel High Viscosity Ratio Applicator with 10 cm Mixing Tip
ABS-10051	ViscoSpray Low Viscosity Ratio Applicator with 3 cm Mixing/Spray Tip
ABS-20000	Fenestrated Delivery Needle
ABS-21000	Tuohy Delivery Needle, 17 Gauge 16"
AR-5995-ABS01	Cart for Centrifuge

An Anticoagulant can be purchased on request.

This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's Directions For Use.

© 2012, Arthrex Medizinische Instrumente GmbH. All rights reserved.

LABS2-003-EN_I

References

1. Borzini P, Mazzucco L: Tissue Regeneration and in Loco Administration of Platelet Derivates: Clinical Outcomes, Heterogeneous Products, and Heterogeneity of Effector Mechanisms. *Transfusion*. 2005; 45: 1759-1767.
2. Edwards D, et al: Transforming Growth Factor Beta Modulates the Expression of Collagenase and Metalloproteinase Inhibitor. *The EMBO Journal*. 1987; 6(7): 1899-1904.
3. Lynch S, et al: Role of Platelet-derived Growth Factor in Wound Healing: Synergistic Effects with other Growth Factors. *Proc. Natl. Acad. Sci. USA*. 1987; 84: 7696-7700.
4. Graziani F, et al: The In Vitro Effect of Different PRP Concentrations on Osteoblasts and Fibroblasts. *Clin Oral Implants Res*. 2006; 17(2): 212-219.
5. Sundman E, Cole B, Fortier L: Growth Factor and Catabolic Cytokine Concentrations Are Influenced by the Cellular Composition of Platelet-Rich Plasma. *American Journal of Sports Medicine*. 2011; 39(10): 2135-2140
6. Kisiday J, et al: Effects of Platelet-Rich Plasma Composition on Anabolic and Catabolic Activities in Equine Cartilage and Meniscal Explants. *Cartilage*. 2012; 3: 245-254
7. Mazzocca A, et al: Platelet-rich plasma differs according to preparation method and human variability. *Journal of Bone & Joint Surgery*. 2012; 94(4): 308-316
8. Mazzocca A, et al: The positive effects of different platelet-rich plasma methods on human muscle, bone, and tendon cells. *The American Journal of Sports Medicine*. 2012; 40(8): 1742-1749
9. Milano G, et al: Repeated Platelet Concentrate Injections Enhance Reparative Response of Microfractures in the Treatment of Chondral Defects of the Knee: An Experimental Study in an Animal Model. *Arthroscopy*. 2012; 28(5): 688-701
10. Cerza F, et al: Comparison Between Hyaluronic Acid and Platelet-Rich Plasma, Intra-articular Infiltration in the Treatment of Gonarthrosis. *The American Journal of Sports Medicine*, online published on October 25, 2012
11. Deans VM, Miller A, Ramos J: A Prospective Series of Patients with Chronic Achilles Tendinopathy Treated with Autologous-conditioned Plasma Injections Combined with Exercise and Therapeutic Ultrasonography. *The Journal of Foot and Ankle Surgery*. 2012; 51(6): 706-710